

Anti-Hk1 Monoclonal Antibody For Wb, If/Icc, Elisa - P19367

Item Number: CM0006712



Introduction

Rabbit monoclonal antibody against human HK1, validated for WB, IF/ICC, and ELISA. Suitable for human, mouse, and rat samples. HK1 catalyzes the first step of glycolysis and plays a role in innate immunity.

[Learn More](#)

Parameter	Specification
Product Name	HK1 Rabbit mAb
Remarks/Alias	HK; HKD; HKI; HXK1; NMSR; RP79; HMSNR; HK1-ta; HK1-tb; HK1-tc; NEDVIBA; hexokinase
Species	Human
GeneID (Human)	3098
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-219 of human HK1 (NP_000179.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19367
Protein Weight	102kDa
Shipping	Ice bag