


Anti-Hk3 Polyclonal Antibody For Wb, If/Icc, Elisa - P52790

Item Number: CM0026389



Introduction

Rabbit polyclonal antibody against human Hexokinase-3 (HK3), suitable for WB, IF/ICC, and ELISA. Detects HK3 in human, mouse, and rat samples, with a predicted band at 99 kDa. Ideal for glycolysis and cancer metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	HK3 Rabbit pAb
Remarks/Alias	HXK3, HKIII, HK3
Species	Human
Gene ID (Human)	3101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human HK3 (NP_002106.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52790
Protein Weight	99kDa
Shipping	Ice bag