

Anti-Hkdc1 Rabbit Polyclonal Antibody For Wb And Elisa - Q2Tb90

Item Number: CM0021288

Introduction

Anti-HKDC1 Rabbit Polyclonal Antibody validated for Western blot and ELISA. Detects human HKDC1, a hexokinase involved in glucose homeostasis and hepatic lipid accumulation. Ideal for metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	HKDC1 Rabbit pAb
Remarks/Alias	RP92; HKDC1
Species	Human
GenelD (Human)	80201
GenelD	80201
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q2TB90
Protein Weight	103kDa
Shipping	Ice bag