

Anti-Hla-B Rabbit Polyclonal Antibody For Wb - P01889

Item Number: CM0028365



Introduction

Rabbit polyclonal antibody against human HLA-B, validated for Western blot. Targets a major histocompatibility complex class I molecule involved in antigen presentation to CD8+ T cells and immune surveillance.

[Learn More](#)

Parameter	Value
Product Name	MITD Rabbit pAb
Remarks/Alias	HLA-B
Species	Human
Gene ID	3106
Immunogen	Synthetic peptide corresponding to amino acids 263-362 of human HLA-B
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P01889
Protein Weight	40KDa
Shipping	Ice bag