

Anti-Hla B7 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P01889

Item Number: CM0006269



Introduction

High-specificity rabbit monoclonal antibody against human HLA B7, validated for Western blot, immunofluorescence, and ELISA applications. Ideal for immune response and MHC I studies. SWISS: P01889.

[Learn More](#)

Parameter	Value
Product Name	HLA B7 Rabbit mAb
Remarks/Alias	AS; HLAB; B-4901; HLA B7
Species	Human
GeneID	3106
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HLA B7 (NP_005505.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P01889
Protein Weight	40kDa
Shipping	Ice bag