


Anti-Hla Class 1 Abc Polyclonal Antibody For Wb, Ihc-P, Elisa - P04439 / P01889 / P10321

Item Number: CM0027809



Introduction

Rabbit polyclonal antibody targeting human HLA-A, HLA-B, and HLA-C. Suitable for Western blot, IHC-P, and ELISA. Ideal for immunological research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	HLA Class 1 ABC Rabbit pAb
Remarks/Alias	AS; HLAB; B-4901
Species	Human
GeneID (Human)	3106
GeneID	3105/3106/3107
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human HLA-A+HLA-B+HLA-C (NP_001229687.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P04439/P01889/P10321
Protein Weight	40 kDa
Shipping	Ice bag