

Anti-Hla-Dmb Rabbit Monoclonal Antibody For Wb, Elisa - P28068

Item Number: CM0003917



Introduction

Rabbit monoclonal antibody against human HLA-DMB (MHC class II antigen DM beta chain). Validated for WB and ELISA. Cross-reacts with human. Molecular weight 29 kDa. Ideal for immunology and antigen presentation research.

[Learn More](#)

Parameter	Value
Product Name	HLA-DMB Rabbit mAb
Remarks/Alias	RING7; D6S221E; HLA-DMB
Species	Human
Gene ID (Human)	3109
Gene ID	3109
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HLA-DMB (P28068).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P28068
Protein Weight	29kDa
Shipping	Ice bag