

Anti-Hla-Dpb1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P04440

Item Number: CM0011604



Introduction

Rabbit monoclonal antibody against human HLA-DPB1. Validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Detects a 29 kDa band. Ideal for MHC class II antigen presentation research.

[Learn More](#)

Parameter	Value
Product Name	HLA-DPB1 Rabbit mAb
Remarks/Alias	DPB1; HLA-DP; HLA-DPB; HLA-DP1B
Species	Human
GeneID (Human)	3115
GeneID	3115
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human HLA-DPB1 (NP_002112.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04440
Protein Weight	29kDa
Shipping	Ice bag