

Anti-Hla-Dqa1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P01909

Item Number: CM0002322

Introduction

Recombinant rabbit monoclonal antibody against human HLA-DQA1 (MHC class II DQ alpha chain). Validated for WB, IHC-P, ELISA, cross-reacts with mouse and rat. Ideal for immunology and antigen presentation research.

[Learn More](#)

Parameter	Value
Product Name	HLA-DQA1 Rabbit mAb
Remarks/Alias	DQA1; DQ-A1; CELIAC1; HLA-DQA; HLA-DQB1; HLA-DQA1*; HLA-DQA1
Species	Human
GeneID (Human)	3117
GeneID	3117
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 155-254 of human HLA-DQA1 (P01909)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01909
Protein Weight	28kDa
Shipping	Ice bag