

Anti-Hla-Dqa1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P01909

Item Number: CM0023875



Introduction

HLA-DQA1 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human and mouse. Targets MHC class II alpha chain involved in antigen processing and presentation. Ideal for immunology and autoimmune disease research.

[Learn More](#)

Parameter	Value
Product Name	HLA-DQA1 Rabbit pAb
Remarks/Alias	DQA1; DQ-A1; CELIAC1; HLA-DQA; HLA-DQB1; HLA-DQA1*; HLA-DQA1
Species	Human
GeneID (Human)	3117
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-216 of human HLA-DQA1 (NP_002113.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01909
Protein Weight	28kDa
Shipping	Ice bag