

Anti-Hla-Dqb1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P01920

Item Number: CM0021248



Introduction

Anti-HLA-DQB1 rabbit polyclonal antibody for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Targets human HLA class II DQ beta 1 chain (MHC II); ideal for immunology and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	HLA-DQB1 Rabbit pAb
Remarks/Alias	IDDM1; CELIAC1; HLA-DQB; HLA-DQB1
Species	Human
GeneID (human)	3119
GeneID	3119
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 33-127 of human HLA-DQB1 (NP_002114.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01920
Protein Weight	30kDa
Shipping	Ice bag