

Anti-Hla-Dr Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P01903

Item Number: CM0006108



Introduction

Rabbit monoclonal antibody targeting HLA-DR alpha chain (P01903). Validated for WB, IF, IHC-P, and ELISA in human samples; ideal for studying antigen presentation, T-cell responses, and immunology research.

[Learn More](#)

Parameter	Value
Product Name	HLA-DR Rabbit PolymAb®
Remarks/Alias	HLA-DRA1
Species	Human
GeneID (Human)	3122
GeneID	3122
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human HLA-DRA (P01903)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag