

Anti-Hla-Dr Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P01903

Item Number: CM0008417



Introduction

Recombinant rabbit monoclonal antibody against human HLA-DRA, validated for WB, IF/ICC, ELISA. Targets MHC class II antigen DR alpha chain, 29kDa, suitable for immunology and immuno-oncology research.

[Learn More](#)

Parameter	Value
Product Name	HLA-DR Rabbit mAb
Remarks/Alias	HLA-DRA1
Species	Human
Gene ID (Human)	3122
Gene ID	3122
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-126 of human HLA-DRA (NP_061984.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag