

Anti-Hla-Dra Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P01903

Item Number: CM0002475



Introduction

High-specificity rabbit monoclonal antibody targeting human HLA-DRA. Validated for WB, IHC-P, and ELISA, with no cross-reactivity to other species. Detects a 29 kDa band. Ideal for immunology and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	HLA-DRA Rabbit mAb
Remarks/Alias	HLA-DRA1; HLA-DRA
Species	Human
GeneID (Human)	3122
GeneID	3122
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HLA-DRA (P01903).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag