

Anti-Hla-Dra Rabbit Polyclonal Antibody For Wb, Elisa - P01903

Item Number: CM0018613



Introduction

Anti-HLA-DRA Rabbit Polyclonal Antibody (P01903) for Western blot and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying MHC class II antigen presentation and T cell biology.

[Learn More](#)

Parameter	Value
Product Name	HLA-DRA Rabbit pAb
Remarks/Alias	HLA-DRA1; HLA-DRA
Species	Human
GeneID	3122
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-216 of human HLA-DRA (NP_061984.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag