

Anti-Hla-Dra Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P01903

Item Number: CM0016137



Introduction

Rabbit polyclonal antibody against HLA-DRA, validated for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Cross-reacts with human and rat. Suitable for immunology and antigen presentation research.

[Learn More](#)

Parameter	Value
Product Name	HLA-DRA Rabbit pAb
Remarks/Alias	HLA-DRA1; HLA-DRA
Species	Human
GeneID (human)	3122
GeneID	3122
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-216 of human HLA-DRA (NP_061984.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag