

## Anti-Hla-Drb3 Polyclonal Antibody For Wb And Elisa - P79483

Item Number: CM0013139



### Introduction

Rabbit polyclonal antibody against human HLA-DRB3 (MHC class II DR beta 3 chain). Validated for WB and ELISA, cross-reacts with Human, Mouse, Rat. Ideal for studying adaptive immunity and antigen presentation.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | HLA-DRB3 Rabbit pAb                                                                                                   |
| Remarks/Alias    | DRB3; HLA-DPB1; HLA-DR1B; HLA-DR3B; HLA-DRB3*; HLA-DRB3                                                               |
| Species          | Human                                                                                                                 |
| Gene ID (Human)  | 3125                                                                                                                  |
| Gene ID          | 3125                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 30-227 of human HLA-DRB3 (NP_072049.2). |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | P79483                                                                                                                |
| Protein Weight   | 30kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |