

Anti-Hla-Drb4 Rabbit Monoclonal Antibody For Wb And Elisa - P13762

Item Number: CM0002462

Introduction

Rabbit monoclonal antibody against human HLA-DRB4 (MHC class II DR beta 4 chain). Validated for WB and ELISA. Presents processed antigens to CD4+ T cells, essential for adaptive immunity studies and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	HLA-DRB4 Rabbit mAb
Remarks/Alias	DR4; DRB4; HLA-DRB; HLA-DR4B; HLA-DRB4
Species	Human
GeneID (Human)	3126
GeneID	3126
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human HLA-DRB4 (P13762).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13762
Protein Weight	30kDa
Shipping	Ice bag