

Anti-Hla-E Polyclonal Antibody For Wb And Elisa - P13747

Item Number: CM0026397



Introduction

Rabbit polyclonal antibody targeting human HLA-E (MHC class I antigen E) for Western blot and ELISA. Validated in human samples, this antibody aids immune surveillance and NK cell studies.

[Learn More](#)

Parameter	Value
Product Name	HLA-E Rabbit pAb
Remarks/Alias	QA1; HLA-6.2; HLA-E
Species	Human
GeneID (Human)	3133
GeneID	3133
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human HLA-E (NP_005507.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13747
Protein Weight	40kDa
Shipping	Ice bag