

Anti-Hla-G Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P17693

Item Number: CM0007203



Introduction

Rabbit monoclonal antibody against human HLA-G (MHC-G). Suitable for western blot, IF/ICC, and ELISA. Ideal for studying immune regulation at the maternal-fetal interface. Cross-reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	HLA-G Rabbit mAb
Remarks/Alias	MHC-G; HLA-G
Species	Human
GeneID (Human)	3135
GeneID	3135
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human HLA-G (NP_002118.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P17693
Protein Weight	38kDa
Shipping	Ice bag