


Anti-Hltf Rabbit Monoclonal Antibody For Wb, Elisa - Q14527

Item Number: CM0012335



Introduction

Rabbit monoclonal antibody against HLTf (SWISS Q14527) for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Targets DNA-dependent ATPase/E3 ligase involved in DNA damage tolerance and chromatin regulation.

[Learn More](#)

Parameter	Value
Product Name	HLTf Rabbit mAb
Remarks/Alias	ZBU1; HLTf1; RNF80; HIP116; SNF2L3; HIP116A; SMARCA3; HLTf
Species	Human
Gene ID (Human)	6596
Gene ID	6596
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 900-1009 of human HLTf (Q14527).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14527
Protein Weight	114kDa
Shipping	Ice bag