

Anti-Hlx Rabbit Polyclonal Antibody For Wb, Elisa - Q14774

Item Number: CM0030233

Introduction

Anti-HLX rabbit polyclonal antibody validated for WB and ELISA. Targets human H2.0-like homeobox protein involved in Th1 cell differentiation. Suitable for embryogenesis and hematopoiesis research.

[Learn More](#)

Parameter	Value
Product Name	HLX Rabbit pAb
Remarks/Alias	HB24; HLX1; HLX
Species	Human
GeneID	3142
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14774
Protein Weight	51kDa
Shipping	Ice bag