

Anti-Hmbs Polyclonal Antibody For Wb, Elisa - P08397

Item Number: CM0024539



Introduction

Rabbit polyclonal antibody targeting human HMBS (porphobilinogen deaminase), suitable for WB and ELISA. Cross-reacts with human and mouse. Essential for heme biosynthesis and porphyrin metabolism research.

[Learn More](#)

Parameter	Value
Product Name	HMBS Rabbit pAb
Remarks/Alias	UPS; PBGD; PORC; PBG-D; HMBS
Species	Human
GeneID (human)	3145
GeneID	3145
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-361 of human HMBS (NP_000181.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08397
Protein Weight	39kDa
Shipping	Ice bag