

Anti-Hmbs Rabbit Monoclonal Antibody For Wb, Elisa - P08397

Item Number: CM0002090



Introduction

HMBS Rabbit monoclonal antibody, validated for WB and ELISA. Cross-reacts with human, mouse, rat HMBS (porphobilinogen deaminase), a heme biosynthesis enzyme (UniProt P08397).

[Learn More](#)

Parameter	Value
Product Name	HMBS Rabbit mAb
Remarks/Alias	UPS; PBGD; PORC; PBG-D; HMBS
Species	Human
Gene ID	3145
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human HMBS (P08397).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08397
Protein Weight	39kDa
Shipping	Ice bag