

Anti-Hmg20A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Np66

Item Number: CM0028630

Introduction

Rabbit polyclonal antibody targeting HMG20A, suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for chromatin regulation and neuronal differentiation studies. - Q9NP66

[Learn More](#)

Parameter	Value
Product Name	HMG20A Rabbit pAb
Remarks/Alias	HMGX1; HMGXB1; HMG20A
Species	Human
Gene ID (Human)	10363
Gene ID	10363
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-347 of human HMG20A (NP_060670.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NP66
Protein Weight	40kDa
Shipping	Ice bag