

## Anti-Hmg20B Rabbit Polyclonal Antibody For Wb, Elisa - Q9P0W2

Item Number: CM0020339

### Introduction

Rabbit polyclonal antibody targeting HMG20B (BRAF35, SOXL), suitable for Western blot and ELISA applications with cross-reactivity to human and mouse samples. Ideal for chromatin, cell cycle, and transcription research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | HMG20B Rabbit pAb                                                                                                 |
| Remarks/Alias    | SOXL; HMGX2; BRAF25; BRAF35; HMGXB2; PP7706; pp8857; SMARCE1r; HMG20B                                             |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 10362                                                                                                             |
| Gene ID          | 10362                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-317 of human HMG20B (NP_006330.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | Q9P0W2                                                                                                            |
| Protein Weight   | 36kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |