

Anti-Hmg20B Rabbit Polyclonal Antibody [Ko Validated] For Wb And Elisa - Q9P0W2

Item Number: CM0021752

Introduction

Rabbit polyclonal antibody targeting HMG20B (BRAF35, SMARCE1r). Validated for WB and ELISA; reacts with human, mouse, and rat. Ideal for cell cycle and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HMG20B Rabbit pAb
Remarks/Alias	SOXL; HMGX2; BRAF25; BRAF35; HMGXB2; PP7706; pp8857; SMARCE1r; 0B
Species	Human
GeneID	10362
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P0W2
Protein Weight	36kDa
Shipping	Ice bag