

Anti-Hmga1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa, Chip-Seq - P17096

Item Number: CM0010366



Introduction

Rabbit monoclonal antibody targeting human HMGA1 (P17096), suitable for WB, IHC-P, IF/ICC, ELISA, and ChIP-seq. Cross-reacts with human, mouse, rat. HMGA1 is a DNA-binding protein involved in chromatin structure and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	HMGA1 Rabbit mAb
Remarks/Alias	HMG-R; HMGY; HMGA1A; HMGA1
Species	Human
GeneID	3159
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-107 of human HMGA1 (P17096).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP-seq
Cross Reactivity	Human, Mouse, Rat
SWISS	P17096
Protein Weight	12kDa
Shipping	Ice bag