

Anti-Hmga2 Rabbit Mab For Wb, Ihc-P, If/Icc, Elisa - P52926

Item Number: CM0007609

Introduction

Rabbit monoclonal antibody recognizing HMGA2 (high mobility group AT-hook protein 2). Validated for WB, IHC-P, IF/ICC, ELISA in human and mouse. Ideal for studying transcriptional regulation and cell cycle.

[Learn More](#)

Parameter	Value
Product Name	HMGA2 Rabbit mAb
Remarks/Alias	BABL; LIPO; SR55; HMGIC; HMGI-C; STQTL9; HMGA2
Species	Human
GeneID (human)	8091
GeneID	8091
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HMGA2 (NP_003474.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52926
Protein Weight	12 kDa
Shipping	Ice bag