

Anti-Hmga2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P52926

Item Number: CM0028977

Introduction

High-quality HMGA2 Rabbit Polyclonal Antibody (P52926), validated for Western blot, immunohistochemistry, and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for transcriptional regulation and chromatin structure research.

[Learn More](#)

Parameter	Specification
Product Name	HMGA2 Rabbit pAb
Remarks/Alias	BABL; LIPO; SRS5; HMGIC; HMGI-C; STQTL9; HMGA2
Species	Human
Gene ID (Human)	8091
Gene ID	8091
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-109 of human HMGA2 (NP_003474.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52926
Protein Weight	12kDa
Shipping	Ice bag