


Anti-Hmgb1 Rabbit Monoclonal Antibody (Ko Validated) For Wb, Ihc-P, If/Icc, Elisa - P09429

Item Number: CM0003189



Introduction

KO validated rabbit monoclonal antibody against human HMGB1 (P09429). Suitable for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Relevant for chromatin biology, inflammation, and autophagy research.

[Learn More](#)

Product Name	[KO Validated] HMGB1 Rabbit mAb
Remarks/Alias	HMG1; HMG3; HMG-1; SBP-1; B1
Species	Human
Gene ID	3146
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HMGB1 (P09429).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09429
Protein Weight	25kDa
Shipping	Ice bag