

Anti-Hmgb2 Polyclonal Antibody For Wb, If/Icc, Elisa - P26583

Item Number: CM0021234

Introduction

Rabbit polyclonal antibody to HMGB2, validated for WB, IF/ICC, ELISA. Detects human, mouse, and rat HMGB2. Ideal for chromatin and transcription research.

[Learn More](#)

Parameter	Value
Product Name	HMGB2 Rabbit pAb
Remarks/Alias	HMG2; HMGB2
Species	Human
Gene ID (Human)	3148
Gene ID	3148
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26583
Protein Weight	24kDa
Shipping	Ice bag