

Anti-Hmgb2 Rabbit Monoclonal Antibody For Wb - P26583

Item Number: CM0011750

Introduction

Rabbit monoclonal antibody targeting human HMGB2 (HMG2), validated for WB, IHC-P, IF/ICC, IP, ELISA, ChIP, and CUT&Tag. Detects HMGB2 in human, mouse, and rat samples. Ideal for chromatin and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	HMGB2 Rabbit mAb
Remarks/Alias	HMG2; HMGB2
Species	Human
Gene ID (Human)	3148
Gene ID	3148
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HMGB2 (P26583).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP, CUT&Tag
Cross Reactivity	Human, Mouse, Rat
SWISS	P26583
Protein Weight	24 kDa
Shipping	Ice bag