

Anti-Hmgb2 Rabbit Polyclonal Antibody For Wb, Elisa - P26583

Item Number: CM0030119

Introduction

Rabbit polyclonal antibody against HMGB2, validated for Western blotting and ELISA. Reacts with human, mouse, and rat HMGB2 (24 kDa). Suitable for chromatin structure, transcription, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	HMGB2 Rabbit pAb
Remarks/Alias	HMGB2; HMGB2
Species	Human
Gene ID (Human)	3148
Gene ID	3148
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human HMGB2 (P26583)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26583
Protein Weight	24kDa
Shipping	Ice bag