


Anti-Hmgb3 Rabbit Polyclonal Antibody For Wb And Elisa - O15347

Item Number: CM0015840



Introduction

High-quality rabbit polyclonal antibody targeting human HMGB3 (high mobility group protein B3, HMG4/HMG-2a). Suitable for Western blot and ELISA; cross-reacts with rat. Ideal for chromatin, transcription, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	HMGB3 Rabbit pAb
Remarks/Alias	HMG4; HMG-4; HMG2A; HMG-2a; HMGB3
Species	Human
GeneID (human)	3149
GeneID	3149
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	O15347
Protein Weight	23kDa
Shipping	Ice bag