

Anti-Hmgcl Polyclonal Antibody For Wb, Ihc-P, Elisa - P35914

Item Number: CM0027585



Introduction

Rabbit polyclonal antibody against human HMGCL, a mitochondrial lyase critical for ketogenesis and leucine catabolism. Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for metabolic disorder research.

[Learn More](#)

Parameter	Value
Product Name	HMGCL Rabbit pAb
Remarks/Alias	HL, HMGCL
Species	Human
GeneID (Human)	3155
GeneID	3155
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human HMGCL (NP_000182.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35914
Protein Weight	34kDa
Shipping	Ice bag