


Anti-Hmgcs1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q01581

Item Number: CM0029658



Introduction

Rabbit polyclonal antibody targeting human HMGCS1, validated in WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for cholesterol metabolism and lipid biosynthesis studies.

[Learn More](#)

Parameter	Specification
Product Name	HMGCS1 Rabbit pAb
Remarks/Alias	HMGCS; HMGCS1
Species	Human
GeneID (human)	3157
GeneID	3157
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 241-520 of human HMGCS1 (NP_002121.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01581
Protein Weight	57kDa
Shipping	Ice bag