


Anti-Hmgcs2 Rabbit Monoclonal Antibody For Wb, If, Ihc-P, Elisa - P54868

Item Number: CM0006959



Introduction

Rabbit monoclonal antibody against human HMGCS2, validated for WB, IF, IHC-P, and ELISA. Cross-reacts with mouse and rat. Useful for ketogenesis and mitochondrial metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	HMGCS2 Rabbit mAb
Remarks/Alias	HMGCS2; Hydroxymethylglutaryl-CoA synthase; mitochondrial; 3-hydroxy-3-methylglutaryl coenzyme A synthase; 3-hydroxy-3-methylglutaryl-CoA synthase 2
Species	Human
Gene ID	3158
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human HMGCS2. (NP_005509.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54868
Protein Weight	50KDa/52KDa/56KDa
Shipping	Ice bag