


Anti-Hmgcs2 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P54868

Item Number: CM0013903



Introduction

Rabbit polyclonal antibody targeting human HMGCs2, validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for ketogenesis, lipid metabolism, and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	HMGCs2 Rabbit pAb
Remarks/Alias	HMGCs2
Species	Human
GeneID	3158
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 269-508 of human HMGCs2 (NP_005509.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54868
Protein Weight	57kDa
Shipping	Ice bag