

Anti-Hmgn1 Polyclonal Antibody For Wb, If/Icc, Elisa - P05114

Item Number: CM0025140

Introduction

Rabbit polyclonal antibody targeting human HMGN1 (HMG14), a non-histone chromosomal protein. Validated for WB, IF/ICC, and ELISA; reacts with human and mouse. Suitable for chromatin and gene expression research.

[Learn More](#)

Parameter	Value
Product Name	HMGN1 Rabbit pAb
Remarks/Alias	HMG14; HMGN1
Species	Human
GeneID	3150
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HMGN1 (NP_004956.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05114
Protein Weight	11kDa
Shipping	Ice bag