

Anti-Hmox1 (Ho-1) Rabbit Monoclonal Antibody For Intracellular Flow Cytometry - P09601

Item Number: CM0008261



Introduction

-conjugated rabbit monoclonal antibody targeting human HMOX1 (HO-1) for intracellular flow cytometry. Validated for apoptosis, oxidative stress, and heme degradation studies. SWISS: P09601.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human HMOX1 mAb
Remarks/Alias	HO-1; HSP32; HMOX1D; bK286B10
Species	Human
GeneID (Human)	3162
GeneID	3162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human Heme Oxygenase 1 (HO-1/HMOX1) (NP_002124.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P09601
Protein Weight	33kDa
Shipping	Ice bag