


Anti-Hmox2 Rabbit Polyclonal Antibody For Wb, Elisa - P30519

Item Number: CM0027903



Introduction

Rabbit polyclonal antibody targeting HMOX2 (heme oxygenase-2) validated for Western blot and ELISA. Reacts with human, mouse, and rat HMOX2. Ideal for heme catabolism and CO signaling research.

[Learn More](#)

Parameter	Value
Product Name	HMOX2 Rabbit pAb
Remarks/Alias	HO-2; HMOX2
Species	Human
GeneID (Human)	3163
GeneID	3163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human HMOX2 (NP_002125.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30519
Protein Weight	36kDa
Shipping	Ice bag