

Anti-Hn1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Uk76

Item Number: CM0006900

Introduction

Rabbit monoclonal antibody against human HN1 (JPT1), suitable for Western blot and ELISA. Recognizes a 16 kDa protein involved in cell cycle regulation, adhesion, and androgen receptor signaling.

[Learn More](#)

Parameter	Value
Product Name	HN1 Rabbit mAb
Remarks/Alias	HN1; ARM2; HN1A
Species	Human
GeneID (human)	51155
GeneID	51155
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-154 of human HN1 (NP_057269.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UK76
Protein Weight	16kDa
Shipping	Ice bag