

Anti-Hnf-4-Alpha Monoclonal Antibody For Wb, Ihc-P, Elisa, And Chip - P49698

Item Number: CM0002164

Introduction

Rabbit monoclonal antibody targeting HNF-4-alpha, a nuclear receptor critical for liver development and circadian rhythm. Validated for WB, IHC-P, ELISA, and ChIP with cross-reactivity in human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	HNF-4-alpha Rabbit mAb
Remarks/Alias	Hnf4; HNF-4; MODY1; Nr2a1; Tcf14; TCF-14; Hnf4alpha; HNF-4-alpha
Species	Mouse
Gene ID (Human)	3172
Gene ID	15378
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P49698
Protein Weight	53kDa
Shipping	Ice bag