


Anti-Hnf1A Rabbit Polyclonal Antibody For Wb, Elisa - P20823

Item Number: CM0029124



Introduction

Rabbit polyclonal antibody against human HNF1α (Hepatocyte nuclear factor 1-alpha), suitable for WB and ELISA. Detects endogenous HNF1α in human, mouse, and rat samples. For research on liver-specific transcription and diabetes.

[Learn More](#)

Parameter	Specification
Product Name	HNF1α Rabbit pAb
Remarks/Alias	HNF1; LFB1; TCF1; HNF4A; MODY3; TCF-1; HNF-1A; IDDM20; HNF1alpha; HNF-1-alpha; HNF1α
Species	Human
GeneID (Human)	6927
GeneID	6927
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-130 of human HNF1α (NP_000536.6).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20823
Protein Weight	67kDa
Shipping	Ice bag