

Anti-Hnf3B/Foxa2 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q9Y261

Item Number: CM0003068



Introduction

Rabbit monoclonal antibody targeting human HNF3 β /FOXA2 (Q9Y261). Suitable for WB, IF/ICC, IP, and ELISA applications. Ideal for research on development, glucose homeostasis, and liver function.

[Learn More](#)

Parameter	Value
Product Name	HNF3 β /FOXA2 Rabbit mAb
Remarks/Alias	HNF3B; TCF3B; HNF-3-beta; HNF3 β /FOXA2
Species	Human
GeneID (human)	3170
GeneID	3170
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HNF3 β /FOXA2 (Q9Y261)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q9Y261
Protein Weight	48kDa
Shipping	Ice bag