

Anti-Hnf4A Polyclonal Antibody For Wb And Elisa - P41235

Item Number: CM0019741



Introduction

Rabbit polyclonal antibody against human HNF4A, validated for WB and ELISA. Suitable for studying transcriptional regulation of hepatic genes, circadian rhythm, and metabolic diseases.

[Learn More](#)

Parameter	Value
Product Name	HNF4A Rabbit pAb
Remarks/Alias	TCF; HNF4; MODY; FRTS4; MODY1; NR2A1; TCF14; HNF4a7; HNF4a8; HNF4a9; NR2A21; TCF-14; HNF4alpha; HNF4A
Species	Human
GeneID (human)	3172
GeneID	3172
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 226-325 of human HNF4A (NP_849180.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P41235
Protein Weight	53kDa
Shipping	Ice bag