


Anti-Hnf4G Polyclonal Antibody For Wb, If/Icc, Elisa - Q14541

Item Number: CM0023749



Introduction

Rabbit polyclonal antibody targeting human HNF4G (NR2A2), validated for Western blot, IF/ICC, and ELISA; cross-reacts with mouse and rat; detects a ~46 kDa band. Ideal for studying nuclear receptor transcription factor function.

[Learn More](#)

Parameter	Value
Product Name	HNF4G Rabbit pAb
Remarks/Alias	NR2A2; NR2A3; HNF4G
Species	Human
GeneID (Human)	3174
GeneID	3174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 366-445 of human HNF4G (NP_004124.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14541
Protein Weight	46kDa
Shipping	Ice bag