


Anti-Hnf4G Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14541

Item Number: CM0012131



Introduction

Rabbit monoclonal anti-HNF4G antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat HNF4G. Ideal for transcriptional regulation and nuclear receptor research.

[Learn More](#)

Parameter	Value
Product Name	HNF4G Rabbit mAb
Remarks/Alias	NR2A2; NR2A3; HNF4G
Species	Human
Gene ID (Human)	3174
Gene ID	3174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 319-408 of human HNF4G5 (NP_004124.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14541
Protein Weight	46kDa
Shipping	Ice bag