


Anti-Hnmt Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P50135

Item Number: CM0020444



Introduction

Rabbit polyclonal antibody against human HNMT. Validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with mouse. Ideal for histamine metabolism and airway response research.

[Learn More](#)

Parameter	Value
Product Name	HNMT Rabbit pAb
Remarks/Alias	HMT; MRT51; HNMT-S1; HNMT-S2; HNMT
Species	Human
Gene ID (Human)	3176
Gene ID	3176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of human HNMT (NP_008826.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50135
Protein Weight	33kDa
Shipping	Ice bag